

## Open Care System Silver Without Phototherapy Manufactured by Desco Medical India

  
  
Decide with Confidence

  
DELIVERING TO THE WORLD

### Open Care System Silver Without Phototherapy



[www.descomedicalindia.com](http://www.descomedicalindia.com)[www.babycareequipments.com](http://www.babycareequipments.com)[www.descoonline.in](http://www.descoonline.in)

### Specification:-

Open Care System Standard are designed to provide advanced warming capabilities with one touch operations

The digital microprocessor based controller comes with a full spectrum of features.

The Open Care System Silver comes complete with automatic manual and baby (Servo) modes.

Two function control models provide more flexibility while optimizing the thermal environment.

Skin temperature display is larger than set temperature display

Skin temperature display can be converted from °C to °F by the push of a button.

A full set of alarm provides safety and convenience, with both audible and visual indicators in the event of:

High Temp. / Low Temp , Set High / Set Low , Probe Failure , Power Failure



# Open Care System Silver Without Phototherapy (MODEL NO : BCOC 106 )



The system automatically starts the heating from 40% in manual mode and the same can be changed through settings upto 99%.

The High & Low temp. alarm settings can be adjusted according to the patient condition.

The Radiant Warmer features an independent ambient temperature probe which maximizes the heating efficiency across the range of operating temperature.

A big LCD is provided for more specific information & visualization of the heater power in manual and servo modes.

The system automatically converts from servo mode to manual mode in case of probe failure thus providing non stop working of the unit.

The spacious baby bed incorporates easily drop down sided for maximum access to the patient.

**Heater Assembly:** The heater assembly located above the center of the support structure consists of an integral examination light and a 650 Watts heating system for quick warming and even distribution of heat on the baby bed.

The neonate can be supplied with a constant oxygen flow through the oxygen humidifier when connected to an oxygen cylinder

**Electrical Specification**

Power supply : 200V/50Htz

Heater Power : 650Watts

Fuse : 5 Amp

IEC : COMPLIANT

We are Manufacturer, Exporters and Suppliers of Open Care System Silver Without Phototherapy at wholesale competitive prices from India

**MODEL NO:**

**BCOC 106** OPEN CARE SYSTEM SILVER WITHOUT PHOTOTHERAPY