

Physio Hand Exercise Table Supplied by Desco Medical India



Specification:-

Hard wood Table with laminated top is fitted with pulleys, Leather loops, Nylon cord & Weights for fingers, Hand, Wrist & Forearm exercises.

Especially for Metacarpophalangeal / Inetrphalangeal joints. Also provided with Supinator / Pronator & Wrist-Circumduction wheel for wrist exercises.

Technical :

Table : Tubular Steel frame table with laminated top, Size : 60cm X 75cm X 75cm high for Finger & Grip Exercises.

Pulleys : 6 (Six) Nylon pulleys, fitted in Steel frame with Powder Coated Finish.

12 Leather Loops with Nylon cord that passes through the pulleys for finger exercises.

Weight Hangers : Six, Weight Hangers with Pulleys underneath the table to hang weights.

Weights : Five Sets, each set consists of 5 weights of one each of 100, 200, 300, 400 and 500 grams. Weights with a hole to hang have chrome finished.

Physio Hand Exercise Table Supplied by Desco Medical India

Specification:-

For Wrist Exerciser :

Wrist Roller : An aluminium roller with handles for supination- pronation & wrist circumduction is provided for wrist exercises.

Wrist Roller is provided with resistance adjustment knob.

Degree of Rotation can be measured and recorded.

Finish : Metal parts with Powder Coating Finish.

We are Exporters and Suppliers of Physio Hand Exercise Table at wholesale competitive prices from India.

MODEL NO:

PTME 401

Rowing Machine cum Sliding Seat Adjustable Base