

Estd. 1966



DELIVERING TO THE WORLD

Semi Electric OT Table General (MODEL NO : OTSE 105)



Semi Electric OT Table General Manufactured by Desco Medical India



Decide with Confidence

Semi Electric OT Table General



DELIVERING TO THE WORLD



INCLUDED Standard Accessories



www.descomedicalindia.com

www.descoonline.in

www.operatingtablesindia.com

Specification:-

- * Stainless Steel Fitting.
- * Detachable head, leg & Pelvic section
- * Base & Cover made of attractive, east to clean Stainless Steel.
- * Sophisticated mechanics provide smooth, step less articulation of the table Top for precise patient positioning.
- * Center Pillar / In built kidney bridge
- * Five Sectional Radio – Translucent.
- * Height lifting electric with remote control.
- * All other functions are manual.



Semi Electric OT Table General (MODEL NO : OTSE 105)



Technical Specs:-

| | |
|-----------------------------------|--------------|
| Length of Table | 1900 mm |
| Width of Table Top | 520 mm |
| Minimum Height (without Mattress) | Min. 750 mm |
| Maximum Height (without Mattress) | Max. 1000 mm |
| Trendelenburg | 25° |
| Reverse Trendelenburg | 25° |
| Lateral Tilt | ±20° |
| Back Section | ±70° |
| Head Section | ±60° |
| Leg Section | 0°-90° |
| Floor Locking System | Manual |

We are Manufacturer, Exporters and Suppliers of Semi Electric OT Table General at wholesale competitive prices from India.

MODEL NO:

OTSE 105 Semi Electric OT Table General