

Estd. 1966



Sodium Heparin Tube Non Vacuum (MODEL NO : PPTU 901)



DELIVERING TO THE WORLD

Sodium Heparin Tube Non Vacuum Manufactured by Desco Medical India



Decide with Confidence

Sodium Heparin Tube Non Vacuum

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

Routine & Special Biochemical Test from Plasma Action Principal Advantage : The anticoagulant heparin activates antithrombins, thus blocking the coagulation cascade and produce a plasma.

Cap Colour : Green

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently.

We are Manufacturer, Exporters and Suppliers of Sodium Heparin Tube Non Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

Sodium Heparin Tube Non Vacuum (MODEL NO : PPTU 901)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 901	12.4 x 75 MM	2 ML	100 Pc.
PPTU 902	12.4 x 75 MM	3 ML	100 Pc.
PPTU 903	12.4 x 75 MM	4 ML	100 Pc.
PPTU 904	12.4 x 100 MM	5 ML	100 Pc.
PPTU 905	12.4 x 100 MM	6 ML	100 Pc.
PPTU 906	12.4 x 100 MM	7 ML	100 Pc.
PPTU 907	12.4 x 100 MM	8 ML	100 Pc.
PPTU 908	16 x 100 MM	9 ML	100 Pc.
PPTU 909	16 x 100 MM	10 ML	100 Pc.
PPTU 910	16 x 100 MM	11 ML	100 Pc.
PPTU 911	16 x 100 MM	12 ML	100 Pc.