

Estd. 1966



Sodium Heparin Tube Vacuum (MODEL NO : PPTU 881)



DELIVERING TO THE WORLD

Sodium Heparin Tube Vacuum Manufactured by Desco Medical India



Decide with Confidence

Sodium Heparin Tube Vacuum

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in

Specifications:-

Routine & Special Biochemical Test from Plasma Action Principal Advantage : The anticoagulant heparin activates antithrombins, thus blocking the coagulation cascade and produce a plasma.

Cap Colour : Green

Max. Vacuum Degree Tolerance : +10 %

Invert the Tube for 6-8 times gently.

We are Manufacturer, Exporters and Suppliers of Sodium Heparin Tube Vacuum at wholesale competitive prices from India.

Estd. 1966



DELIVERING TO THE WORLD

Sodium Heparin Tube Vacuum (MODEL NO : PPTU 881)



MODEL NO:	OUTER DIA X LENGTH	CAPACITY ML	PACK SIZE
PPTU 881	12.4 x 75 MM	2 ML	100 Pc.
PPTU 882	12.4 x 75 MM	3 ML	100 Pc.
PPTU 883	12.4 x 75 MM	4 ML	100 Pc.
PPTU 884	12.4 x 100 MM	5 ML	100 Pc.
PPTU 885	12.4 x 100 MM	6 ML	100 Pc.
PPTU 886	12.4 x 100 MM	7 ML	100 Pc.
PPTU 887	12.4 x 100 MM	8 ML	100 Pc.
PPTU 888	16 x 100 MM	9 ML	100 Pc.
PPTU 889	16 x 100 MM	10 ML	100 Pc.
PPTU 890	16 x 100 MM	11 ML	100 Pc.
PPTU 891	16 x 100 MM	12 ML	100 Pc.