

Estd. 1966



DELIVERING TO THE WORLD

Syringe Pump Colour Display (MODEL NO : IESY 105)



Syringe Pump Colour Display Supplied by Desco Medical India



Decide with Confidence

Syringe Pump Colour Display

Estd. 1966



DELIVERING TO THE WORLD



www.descomedicalindia.com

www.descoonline.in



Syringe Pump Colour Display (MODEL NO : IESY 105)



Specification:-

A syringe driver or syringe pump is a small infusion pump (some include infuse and withdraw capability), used to gradually administer small amounts of fluid (with or without medication) to a patient or for use in chemical and biomedical research.

Unique Intelligent Self-Recognition technology of syringes

KVO, Bolus and Anti-Bolus

Applicable for 5,10,20,30,50 (60) ml Syringes

Compatible for various brands of syringes

Various Infusion modes to meet different clinical requirements

Various visual and audible alarms

Rotatable pole Clamp

Numeric key design, easy to operate

Large Colorful LCD Screen

Real Time Display, ensure the infusion information to be recorded correctly

3 Adjustable occlusion levels, dynamically display the pressure

Night Mode

Double CPU ensuring safe infusion

We are Exporters and Suppliers of Syringe Pump Colour Display at wholesale competitive prices from India.

MODEL NO

IESY 105

Syringe Pump Colour Display