



# Ventilator Vertical Silver (MODEL NO : APVE 121)



Decide with Confidence

## Ventilator Vertical Silver



DELIVERING TO THE WORLD

[www.descomedicalindia.com](http://www.descomedicalindia.com)

[www.anesthesiaproductsindia.com](http://www.anesthesiaproductsindia.com)

[www.descoonline.in](http://www.descoonline.in)

## Ventilator Vertical Silver Supplied by Desco Medical India

### Specification:-

Ventilator Vertical Silver is a multi-functional Ventilator and is applicable for children and adults in sickroom, ICU or of emergency.

Display Mode : High-definition 5.7" LCD screen display

Ventilation Mode : IPPV, SIPPV, SIMV, MANUAL

Ventilation Function : PEEP, SIGH, IRV

Ventilation Parameters

Tidal volume 50 ~ 1500ml

Rate 2 ~ 99bpm

SIMV rate 2 ~ 20bpm

I:E 2:1 ~ 1:8

Inspiratory trigger pressure -1.0 ~ 2.0 kPa

PEEP 0 ~ 2.0 kPa

Estd. 1966



DELIVERING TO THE WORLD

# Ventilator Vertical Silver (MODEL NO : APVE 121)



- Pressure range 1.0~6.0 kPa
- SIGH 1.5 times the inspiratory time
- Parameters for ventilation monitoring tidal volume, minute ventilation volume, IPPV rate, SIMV rate, total respiratory rate, I/E, peak pressure of airway, pressure - time waveform, flow rate - time waveform, PEEP, inspiratory trigger pressure
- Security Alarm System
  - Airway pressure alarm Upper limit setting range 1.0~6.0 kPa
  - Low limit setting range 0.4~2.0 kPa
  - Per-minute ventilation volume alarm Upper limit setting range 3.0~30L/min.
  - Low limit setting range 1.0~10L/min
  - Sustained high-pressure alarm It will give alarm when stress have consistently been higher than 2.5 kPa
  - Suffocation alarm It will give sound and light alarm if there is no tidal volume input for 15 seconds
  - Power alarm
  - Power AC 220V±10% 50Hz±1Hz UPS, Storage battery
  - Mechanical Arms

We are Exporters and Suppliers of Ventilator Vertical Silver at wholesale competitive prices from India

**Model No:**  
**APVE 121 Ventilator Vertical Silver**